

Before starting

1. Bring samples and buffers to RT.
2. Set thermomixer to 56°C.

Procedure

1. Pipette 1-100 μ l into 1.5ml tube.
2. Make up volume to 100 μ l with ATL buffer.
3. Add 10 μ l proteinase K.
4. Add 100 μ l buffer AL, mix by vortexing.
5. Incubate at 56°C for 10 minutes.
6. Centrifuge briefly.
7. Add 50 μ l ethanol (100%), vortex. Incubate at RT for 3 minutes.
8. Centrifuge briefly.
9. Transfer the entire lysate to Qiagen minielute column which is placed in a 2 ml collection tube.
10. Add 500 μ l of AW1 buffer. Centrifuge at 6000 g for 1 minute. Discard the flowthrough tube.
11. Add 700 μ l of buffer AW2. Centrifuge at 6000 g for 1 minute. Discard the flowthrough tube.
12. Add 700 μ l of ethanol. Centrifuge at 6000 g for 1 minute. Discard the flowthrough tube.
13. Centrifuge at full speed for 3 minutes to dry the membrane.
14. Place column in a fresh 1.5 ml tube (doesn't come with kit). Open the lid and let it dry at RT.
15. Put 20-100 μ l of ATL buffer to the center of the membrane. Incubate for 10 minutes.
16. Centrifuge at full speed for 1 minute.

